

# Quick Start Guide to using the P137/P138 Care Pager

Press for 3 sec to turn on.  
If no response then charge battery for at least 5 minutes



Charger can be left permanently connected or it can boost charge from empty in 8 hours

## Typical screens in normal operation



Shows battery state on bars

Name of alarm will flash up when safety signal received. Bars show strength of signal



Indicates there is an alarm system in the area, not coupled to this pager. This may be OK or an error so check with the supervisor

*The Alert-it system has been designed with due regard to reliability and integrity. While it offers a highly vigilant monitoring method, it is always possible that a distress condition can go undetected for a variety of reasons (including malfunction) and in life threatening situations it is advisable to use the Alert-it system in conjunction with additional monitoring techniques (e.g. video). Neither the manufacturer nor its agent can accept legal responsibility to provide a system that is infallible. The carer is responsible for assessing the risks of using this equipment and any settings pertaining to it.*

**Turn Off**  
Enter MENU (see overleaf)  
Scroll right to Turn Off & Select  
Do not remove batteries as this will damage the connections eventually

### Alarm active



Radio contact with an alarm fitted with "Safelink" security has been lost. This MUST be corrected immediately



Typical alarm indication, with four levels: Urgent, Help, Assist, Fault


1. Attend the user,
2. Reset the alarm on the user's equipment
3. and the pager will return to normal operation automatically

Try walking towards the user's room. If the alarm cancels then there is a problem with the building or distance required. Contact iTs Designs who can offer a signal booster. If the alarm does not cancel then:

1. Check the user equipment is turned on.
2. Check the user equipment is transmitting a regular safety signal (typically the green power light will flash every 8-10 seconds)
3. Contact iTs Designs for further advice

Full adjustment details are found in handbooks available on:  
**[www.alert-it.co.uk/handbooks/](http://www.alert-it.co.uk/handbooks/)**  
Or by phoning Alert-iT

*PTO for setup options*



**Support**

For technical support please fax or EMail:  
 HELP: 0845 2179951  
 FAX : 0845 2179953  
 Support@itsdesigns.co.uk

Designed by:  
 ITs Designs Ltd  
 Leicester  
 LE9 9FE UK

...using technology to care for carers