

Sleep easy.

**Providing users &
carers peace of mind
at night**

**Increases user
confidence and carer
givers awareness.**

Innovative Assistive Technology
For Personal Professional
Care Assistance





Part of the Picomed Group



Class 1 Medical devices – Support for Adults and Children

Product selector

	Companion Mini	Pulse Companion
Fine Movements		
Pulse		

Introducing the Companion Mini range

The Companion Mini is the small, smart, and sensitive unit that can alert the caregiver to an increase in fine movement. It's small, but it's tough and it's completely cable-free. Suitable for children and adults. Designed to detect movement from beneath the user's mattress.

Brilliantly simple to install and use, the Companion Mini's sensor also makes it one of the smallest and most discreet systems yet. There are no cables to connect, and the device is ready in seconds, making it equally useful for use away from home. The unit has a battery life of up to one year, keeping maintenance to a minimum.

Companion Mini For Your Home

A compact, robust device, the Alert-iT Companion Mini detects fine movements from beneath a mattress. Alarm notifications are sent to a Plesio pager using Alert-iT's failsafe radio technology.



Companion Mini Pro The Option for your Healthcare Facility

Designed for use in professional healthcare environments, the Companion Mini Pro is used in conjunction with Alert-iT's advanced pagers (8, 64 or 128 monitor memory). A remote reset unit is supplied for installation in the vicinity of the user, this ensures that the caregiver attends to clear any alarm notifications.



Companion Mini Connect

The Companion Mini Connect is designed to link up with third-party Nurse call and Telecare systems, offer Alert-iT's sensing technology and alarming your preferred call system.



Pulse Companion (MkII)

The Pulse Companion uses a discreet and comfortable wireless sensor, and raises an alarm notification when the detected heart rate reading rises or falls past the user set parameters.

The minimum and maximum sensor readings are recorded throughout the night and with the aid of the Alert-iT App. This helps the caregiver in personalising the alarm parameters for their unique users' requirements. Allows for much less intrusive care, reduces false alarms, and improving the user's independence and dignity.

The caregiver can finely adjust the alarm limits on how fast or slow the users heart rate can go, along with the duration they are allowed to be above or below these set limits. Various safety features are present to detect if the sensor has been removed, locked up in anyway, or simply needs its battery charging.

If the Pulse Companion indicates a possible alerting condition, the caregiver will be notified by either an Alert-iT pager or nurse call/telecare system.

Benefits

- ⦿ Improved independence for the user and reassurance for the caregiver.
- ⦿ Simple to setup via the Alert-iT app.
- ⦿ Peace of mind from Failsafe technology.
- ⦿ 3-year warranty (1 year for wireless sensor).
- ⦿ No license or fees.

Features

- ⦿ Configuration app available on Android and iOS (free).
- ⦿ Fully Personalised to the user.
- ⦿ Battery Powered.



* Smart phone not included.

P163 Bed Alarm

The Alert-iT P163 Bed Alarm is used to support a Person whilst in bed or chair, alarming if they leave.

The Bed Alarm uses a range of sensors to suit a wide range of situations, by placing them on or under the mattress, on a chair or on the floor. The P163 automatically alarms differently if it detects a bed and chair sensor or floor mat.

The P163 Bed alarm has an easy to press button located on the front, so a Person can call for help, whilst also having the capability for the carer to escalate the call to urgent if they require assistance.

The Caregiver is notified of an event via either the Alert-iT failsafe pager or connecting to an existing Nurse Call System.

Having a dual input gives you the flexibility to give a user more ways to call for assistance, for example, input 1 can be timed for bed occupancy while input 2 can have an easy to reach wired call button. Alternatively, if there are two beds or chairs near to each other it can monitor two users at the same time. The P163 is specially designed so it can pick up each Person individually and send separate alarms.



Improving independence and dignity.

The Bed Monitor allows the caregiver to offer less intrusive care, improving the person's independence and dignity.

Sleep Easy Occupancy Sensors

P143 Bed Pads

P143CAB Large bed mat (740x400mm)

A larger than standard sensor mat that helps to reduce false alarms due to its large surface area increasing the contact with the user.

P143CBB Standard bed mat (740x250mm)


A standard size to offer a more cost-effective solution.



P150C Floor Mats

The P150C, is a robust washable very low slip floor mat

Can be connected to most nurse call and Telecare systems with specified lead (sold separately).



P149AD Chair Mat

The P149AD Chair Mat is a large chair mat that can be fitted under, on top or inside an existing cushion.



P144B Bed Weight Sensor

The P144B Bed weight sensor is a unique device that detects someone exiting the bed without contact with the user. Can be used on beds with air low and pressure relieving mattress. Can be used with Telecare systems. 3 year warranty.



P149BAB Chair Ribbon

The P149BAB Chair Ribbon is an under a cushion sensor. Sensor is very flexible and water resistant. 3 year warranty.



P143GAB Bed Ribbon


The P143GAB bed ribbon is an under-mattress sensor. Sensor is very flexible and water resistant. 3 year warranty.



Sleep Easy Moisture Sensors

P142A Cotton Sheet

- Size:** Sensing area size 900x600mm
(total width 1600mm)
- Absorbency:** 5ml before alerting
- Cleaning:** Machine washable recommended
40 degrees (no fabric conditioner)
- Detection:** Enuresis, excessive perspiration,
salivation, and vomit.



P142E Absorbent sheet

- Size:** Sensing area size 900x600mm
(total width 1600mm)
- Absorbency:** 300ml before alerting (500ml in total)
- Cleaning:** Machine washable recommended
40 degrees (no fabric conditioner)
- Detection:** Enuresis.



P142F Pillow Case Sensor

- Size:** Sensing area size 600x450mm
- Absorbency:** 5ml before alerting
- Cleaning:** Machine washable recommended
40 degrees (no fabric conditioner)
- Detection:** Enuresis, excessive perspiration,
salivation, and vomit.






Innovative Assistive Technology For Personal Professional Care Assistance

3-Year Warranty | Top Rated Reliability | ISO Accredited | Over 30 Years Experience

Contact us for a free no obligation chat about
our sleep easy Assistive Technology range

Alert-iT Care Alarms
Fernie House, Unit 3
Jackson Street
Coalville Business Park
Leicester. LE67 3NR
United Kingdom

Social:

-  AlertitCareAlarms
-  Alert_iT_Care
-  Alert-iT

+44 1530 239 900

www.alert-it.co.uk

