



Solutions

TAILORED SOLUTIONS, MADE FOR YOU



**FALLS
MANAGEMENT**



**EPILEPSY
MONITORING**



**DEMENTIA
CARE**



**INCONTINENCE
SUPPORT**



Plesio Care 1 to 1 Solutions



THE GUARDIAN SOLUTION



The Guardian monitor is our most advanced system, making it ideal for demanding situations supporting tonic/clonic seizures and the risk of SUDEP. The Guardian is able to detect a wide range of seizure types and features an integral microphone for detecting vocal seizures. Furthermore, it also offers the potential to monitor seizures indicated by incontinence, vomiting and bed vacuation through a range of accessories. Carers are notified of all alarms via a handheld Alert-iT pager.

Product Code **S1020B**

INCLUDES PLESIO PAGER



THE COMPANION SOLUTION



The Companion Monitor is an effective solution for supporting tonic/clonic seizures as well as a variety of other situations resulting in movement, sound, vomiting or incontinence.

Product Code **S1026C**

Movement and Sound associated seizures.

INCLUDES PLESIO PAGER & MICROPHONE



Product Code **S1026D**

Movement and vomiting or moisture associated seizures.

INCLUDES PLESIO PAGER & 2X COTTON SHEETS



THE COMPANION MINI



The Companion Mini is currently the smallest monitor of its kind available. The digital movement sensor – small enough to fit in the palm of your hand - is our most discreet system yet and very simple to install and use. Offering reliability and peace of mind for supporting tonic-clonic seizures, its compact dimensions and cable-free installation make it ideal for use away from home on weekends away or holidays.

Product Code **S1043A**

INCLUDES PLESIO PAGER





FLOOR MAT SOLUTION



Alert-iT floor mats help you reduce instances of falls by sending an alarm as soon as a person is present on the mat. Slim and unobtrusive, they also feature robust construction and an anti-slip surface for durability with safety. All mats are made from an easy-to-clean, waterproof material to help you with infection control.

Product Code **S1034A**

INCLUDES PLESIO PAGER

Other items that may be required:

Product Code **P159B**

Wall/Bed Bracket Kit for P163.



RADIO PIR MOTION SENSOR



The Alert-iT battery-operated infra-red radio motion sensor detects someone leaving their bed before they become weight bearing. A discreet device, it's controlled by a simple keyfob style-remote control, allowing you to turn the alarm on or off before when entering user's room so they won't be disturbed.

Product Code **S1036B**

INCLUDES PLESIO PAGER



THE CALL BUTTON SOLUTION



The call button can be worn on the wrist or round the neck on a lanyard quick release. It has a large easy press coloured button with audio able and visual reassurance when pressed. Alerts a pager with a range of up to 400 metres.

Product Code **S1032BB**

INCLUDES PLESIO PAGER

The Call/Fall badge is worn around the neck or pinned to the clothing of the upper body. Operates in a similar manner to the call badge, yet has the added benefit of being able to detect that the wearer may have fallen.

Product Code **S1032AB**

INCLUDES PLESIO PAGER





THE BED OCCUPANCY SOLUTION

The Alert-it Bed Occupancy Solution detects when someone leaves the bed and passes an alarm to an Alert-iT handheld pager. As a unique feature; once the alarm has been raised, the alarm will reset automatically if the user returns to the bed within an allotted amount of time, highlighting that they are safe.

Product Code **S1034B**

INCLUDES PLESIO PAGER



Other items that may be required:

Product Code **P159B**

Wall/Bed Bracket Kit for P163.



THE CHAIR OCCUPANCY SOLUTION

The Alert-iT Chair Occupancy Solution detects when someone leaves the chair and passes an alarm to a handheld pager. As a unique feature; once the alarm has been raised, the alarm will reset automatically if the user returns to the chair within an allotted amount of time, highlighting that they are safe.

Product Code **S1034E**

INCLUDES PLESIO PAGER





AUDIBLE BED OCCUPANCY MAT



The Alert-iT Audible Bed Occupancy Mat lets carers know if someone has vacated their bed by raising an audible alarm. The alarm can be set to sound immediately or after a timed interval.

Product Code **S1042B**



AUDIBLE CHAIR OCCUPANCY MAT



The Alert-iT Audible Chair Occupancy Mat lets carers know if someone has vacated their chair by raising an audible alarm. The alarm can be set to sound immediately or after a timed interval.

Product Code **S1042C**



THE DOOR ACTIVITY SOLUTION



The door activity sensor offers many unique features, not only will it alarm the Plesio pager if the door has been opened, but it can be temporarily disarmed while the carer is at the door with automatic re-arming when the door is closed. The unit can only be disabled with a secure key.

Product Code **S1036C**

INCLUDES PLESIO PAGER



Wired Solutions



PIR MOTION DETECTOR

A discreet motion detector that connects directly to your existing nurse call system to raise an alarm as soon as the user's foot descends to the floor. An ideal alternative where floor mats are deemed unsuitable.

PIR powered by 4 AA batteries.

Please note Nurse call leads are included in the RRP, please quote type of system when ordering if required.

Product Code P162AB



FLOOR MAT

A large robust washable low slip floor mat, can be connected to most nurse call systems with specified lead (sold separately).

Product Code P150C



FLOOR MAT CARPET STYLE

A robust carpet floor mat, can be connected to most nurse call systems with specified lead (sold separately).

Product Code P150CB



THE COMPANION



The Companion Monitor is ideal for support of Tonic/Clonic seizures as well as situations that manifest sounds or moisture. It can also warn of possible danger if the user gets out of bed, or stays out too long. Perfect for domestic or professional use, the Companion plugs directly into your existing nurse call system (lead sold separately).

Product Code **S1027AB**

Detects spasm and sound (please note you will require a microphone if detecting sound associated with the seizure).

Product Code **P158A**

Microphone.

Product Code **S1027BB**

Detects spasm and moisture (please note you will cotton sensing sheets to detect moisture associated with the seizure).

Product Code **P142A**

Cotton sensing sheets.



THE GUARDIAN



The Guardian monitor is our most advanced system, making it ideal for demanding situations supporting tonic/clonic seizures and the risk of SUDEP. The Guardian is able to detect a wide range of seizure types and features an integral microphone for detecting vocal seizures. Furthermore, it also offers the potential to monitor seizures indicated by incontinence, vomiting and bed vacation through a range of accessories.

Product Code **S1017A**

Other items that may be required:

Product Code **P142A**

Cotton sensing sheets (to detect moisture).

Product Code **P143C**

Bed mat (to detect bed occupancy).

Alert-iT have distributors throughout Europe who are happy to assist you with our products



Picomed AS

Distributor in Norway (N) and Sweden (SE)

Contact Number: +0047 37119950

Email Address: post@picomed.no

Website: www.picomed.no



Astrid Leisner and Son Ltd

Distributor in Denmark (DK) and Germany (D)

Contact Number: +45 63 71 30 50

Email Address: info@leisner.dk

Website: www.leisner.dk



Tunstall Emergency Response Limited

Distributor in Ireland (IRL)

Contact Number: +353 1850 247 999

Email Address: sales@emergencyresponse.ie

Website: www.tunstallemergencyresponse.ie

If your country is not listed, please contact Alert-iT direct – we'll be happy to help.



Alert-iT Care Alarms

Contact Number: +0044 1530 239900

Email Address: enquiries@alert-it.co.uk

Website: www.alert-it.co.uk

Alert-iT are looking to expand their distributor network,
register your interest at enquiries@alert-it.co.uk

Every solution comes with complete peace of mind...

All Alert-iT technology (excluding bed mats) comes with 3 year warranty and designed for simple installation and use. Our support desk will offer help on the use of our technology if required. If these solutions below don't meet your requirements please call our sales team as we may be able advise and create an individual solution.