Companion Mini

The Companion Mini is the latest addition to our fantastic range of epilepsy monitors and alarms. The digital movement sensor is our most discrete system yet and can be applied in some of the most challenging situations when supporting tonic-clonic seizures.

Currently the smallest monitor of its kind on the market, the device sits comfortably in the palm of your hand. Designed with the highest performance sensors the Companion Mini is able to identify seizure induced movements whilst ignoring sleep activities before alerting a user to irregularities via the Alert-iT Safelink Radio System.

Designed to discreetly monitor an individuals bed-time movements, this user friendly system is simple to set up with little or no adjustment required.

A New Digital Design

Did you know the Companion Mini is the first of its kind?

The fully digital system has been designed to provide a singular solution to support individuals with Epilepsy. With the subtraction of wires associated with other monitors the user is able to freely transport the device with them; whether for nights away or to take on holiday.

Key Features of the Companion Mini:

- Long Lasting Battery The Companion Mini is the first battery powered monitor boasting a life of at least 1 year. This means that there are no protruding wires and further increases the ease of portability.
- Discrete Size The size of the Companion Mini is significantly smaller than any other monitor or alarm.
 Despite its delicate size the monitor continues to provide solid support and peace of mind for individuals. Able to fit in the palm of your hand the monitor is ideal for travelling.
- Casing The Companion mini has been designed with a Biomaster coated plastic case which robust and tamper proof.
- Pager The Companion Mini comes complete with a pager which lasts approximately 3 days but can still be in use when charging.
- A USB cable to link your monitor to a specialised computer programme supplied free of charge is supplied for initial programming.
- Fail safe radio with a 450 metre range
- Ultra-sensitive 1mG movement sensor
- 3 year warranty



How Does It Work?

Easy to install, the Companion Mini provides peace of mind for monitoring the risks and activities for those who have epilepsy.

The monitor is placed under the mattress where it constantly monitors your movements. The Companion Mini is designed to be placed and secured in the centre of the bed between the base and the mattress.

When in place the Companion Mini monitors activity until movements associated with seizures are recognised. Verified and ongoing seizures then result in the transmission of alerts to a hand help pager.



Suitability

Epilepsy is one of the most common neurological conditions in the world affecting approximately 600,000 people alone here in the UK.

Did you know that of those in the UK with epilepsy, a third have idiopathic generalised epilepsy. Of this number approximately 60% of individuals, that's 91,200 people, will experience primary generalised tonic-clonic seizures.

Considering these statistics the Companion Mini has been designed to provide peace of mind for individuals and carers alike in the domestic market.

