



Solutions

TAILORED SOLUTIONS, MADE FOR YOU



Epilepsy



Incontinence



Dementia



Falls





Contents

Section 1:

Plesio Care 1 to 1 Solutions

Section 2:

Wired Solutions

Plesio Care 1 to 1 Solutions



THE GUARDIAN SOLUTION



The Guardian Monitor is our most advanced system and can be applied in some of the most demanding situations when supporting tonic/clonic seizures and SUDEP.

Solution Product Code **S1020B**

£680

INCLUDES PLESIO PAGER



THE COMPANION SOLUTION



The Companion Monitor is an effective solution for supporting tonic/clonic seizures as well as a variety of other situations resulting in movement, sound, vomiting or incontinence.

Solution Product Code **S1026C**

Movement and Sound associated seizures.

£390

INCLUDES PLESIO PAGER & MICROPHONE



Solution Product Code **S1026D**

Movement and vomiting or incontinence associated seizures.

£500

INCLUDES PLESIO PAGER & 2X COTTON SHEETS



THE COMPANION MINI SOLUTION



The Companion Mini, currently the smallest monitor of its kind on the market. The digital movement sensor is our most discrete system yet and can be applied in some of the most challenging situations when supporting tonic-clonic seizures. Very simple to install and operate.

Solution Product Code **S1043A**

£220

INCLUDES PLESIO PAGER





THE CALL BUTTON SOLUTION

The call button can be worn on the wrist or round the neck on a lanyard quick release. It has a large easy press coloured button with audio able and visual reassurance when pressed. Alerts a pager with a range of up to 400 metres.

Solution Product Code **S1032BB**

£165

INCLUDES PLESIO PAGER



The Call/Fall badge is worn around the neck or pinned to the clothing of the upper body. Operates similar to the call badge with the added feature of detecting the wearer may have fallen.

Solution Product Code **S1032AB**

£185

INCLUDES PLESIO PAGER



THE BED OCCUPANCY SOLUTION

A bed occupancy solution detects the person leaving the bed and raises an alarm to the Plesio pager. A unique feature, if the user returns to the bed themselves the alarm will reset automatically highlighting the person is safe. Supplied with a large bed mat to minimize false alarms.

Solution Product Code **S1034B**

£220

INCLUDES PLESIO PAGER



Other items that may be required:

Solution Product Code **P159B** £11.50

Wall/Bed Bracket Kit for P163.



THE CHAIR OCCUPANCY SOLUTION

A chair occupancy solution detects the person leaving the chair and raises an alarm to the Plesio pager. A unique feature, if the user returns to the chair the alarm will reset automatically highlighting the person is safe. Supplied with a large chair mat to minimize false alarms.

Solution Product Code **S1034E**

£243

INCLUDES PLESIO PAGER



Other items that may be required:

Solution Product Code **P159B** £11.50

Wall/Bed Bracket Kit for P163.



THE FLOOR MAT SOLUTION

Our floor mats are robust low-slip mats that do not need over mats for protection or safety. Once the mat is stepped on an alarm is raised to the Plesio pager. The detector unit also has a call button that the user can use as a basic call system.

Solution Product Code **S1034A**

£280

INCLUDES PLESIO PAGER



Other items that may be required:

Solution Product Code **P159B** £11.50

Wall/Bed Bracket Kit for P163.



THE DOOR ACTIVITY SOLUTION

The door activity sensor offers many unique features, not only will it alarm the Plesio pager if the door has been opened, but it can be temporarily disarmed while the carer is at the door with automatic re-arming when the door is closed. The unit can only be disabled with a secure key.

Solution Product Code **S1036C**

£185

INCLUDES PLESIO PAGER



THE DISCREET DETECTION SOLUTION

A discreet detection unit with an infra-red motion detector remote from the bed. The alarm is raised as soon as the user's foot descends to the floor. The same sensor can be positioned to detect them sitting up, leaving a room, or entering a hazardous area. It can also detect door usage, where a door may be left open. It's an alternative to the door sensor and floor mats which can be a trip hazard.

Solution Product Code **S1036B**

£188

INCLUDES PLESIO PAGER



THE COMPLETE SOLUTION

A complete solution which bring together the door sensors and infra-red motion detector. The solution consists of 2x Door, 1x Infra-red motion sensor and our 8 monitor pager which offer a descriptive location and type of device which has been activated. An idea solution to offer independence for the user and peace of mind for the carer.

Solution Product Code **S1036A**

£350

Wired Solutions



WIRED PIR MOTION DETECTOR



A discreet motion detector. The alarm is raised as soon as the user's foot descends to the floor on an Alert-iT pager (radio) or an existing Nurse call system (wired). It's an alternative to the floor mats (if identified as unsuitable).

PIR powered by 4 AA batteries.

Please note Nurse call leads are included in the RRP, please quote type of system when ordering if required.

Solution Product Code P162AB

£88



FLOOR MAT



A large robust washable low slip floor mat, can be connected to most nurse call systems with specified lead (sold separate).

Solution Product Code P150C

£82



FLOOR MAT CARPET STYLE



A robust carpet floor mat, can be connected to most nurse call systems with specified lead (sold separate).

Solution Product Code P150CB

£100



WIRED BED VACATION MONITOR WITH MAT



A Bed occupancy solution detects the person leaving the Bed and raises an alarm to your nurse call system. Supplied with a large bed mat to minimize false alarms.

Radio Monitor powered by 2AA battery's.

Solution Product Code **S1043BB**

£130

Other items that may be required:

Solution Product Code **P159B** **£11.50**

Wall/Bed Bracket Kit for P163.



AUDIBLE BED VACATION MONITOR INCLUDING LARGE BED MAT



A bed occupancy solution detects the person leaving the bed and raises an audible alarm. A unique feature, if the user returns to the chair themselves the alarm will reset automatically highlighting the person is safe. Supplied with a large bed mat to minimize false alarms.

Radio Monitor powered by 2AA battery's.

Solution Product Code **S1042B**

£99

Other items that may be required:

Solution Product Code **P159B** **£11.50**

Wall/Bed Bracket Kit for P163.



THE COMPANION

The Companion Monitor is ideal for support of Tonic/Clonic seizures as well as situation that manifest sounds, vomiting or incontinence. It can also warn of possible danger if the user gets out of bed, or stays out too long. It is perfect for domestic or professional use. Plugs direct in to the existing call system (Lead sold separate).

Solution Product Code S1027AB £225

Detects spasm and sound (please note you will require a microphone if detecting sound associated with the seizure).

Solution Product Code P158A £15

Microphone.

Solution Product Code S1027BB £250

Detects spasm and moisture (please note you will cotton sensing sheets to detect moisture associated with the seizure).

Solution Product Code P142A £55 each

Cotton sensing sheets.



THE GUARDIAN

The Guardian Monitor is our most advanced system and can be applied in some of the most demanding situations when supporting tonic/clonic seizures and SUDEP. Careful design ensures high performance sensors to monitor bed movement, incontinence, vomiting, bed vacation and sounds. The movement detector, which has the capability of supporting those with both tonic/clonic and clonic only seizures, will also continue to monitor shallow bed movements, measuring movements which occur due to breathing and muscle activity the Guardian Monitor will sound an alarm if these movements should cease supporting the risk of SUDEP.

Solution Product Code S1017A £600

Other items that may be required:

Solution Product Code P142A £55 each

Cotton sensing sheets.

Solution Product Code P143C £61

Bed mat (to detect bed occupancy).

Every solution comes with complete peace of mind...

All Alert-iT technology (excluding bed mats) comes with 3 year warranty and designed for simple installation and use. Our support desk will offer help on the use of our technology if required. If these solutions below don't meet your requirements please call our sales team as we may be able advice and create an individual solution.

www.alert-it.co.uk | 0845 217 9952 | enquires@alert-it.co.uk