

## **EU Declaration of Conformity**

Manufacturer:

**European Authorised Representative:** 

iTs Designs Ltd.

(T/A: Alert-iT Care Alarms)

Unit 3 Ground Floor,

Coalville Business Park,

Leicestershire,

LE67 3NR, England.

Tel: +44 (0) 1530 239 900

Leisner ApS

Korsvangcentret

DK-5610 Assens

Denmark

Tel: +45 6371 3050

Single Registration Number (SRN): DK-AR-000027183

Single Registration Number (SRN): GB-MF-000027859

Description:	The P140 is an under mattress movement sensor for use with Alert-iT's bedside monitors.		
Directives:	2017/745	Medical Device Regulation (MDR)	
2	2011/65/EU + 2015/863	Permitted Materials (RoHS)	
	2012/19/EU	Waste Electrical and Electronic Equipment Directive (WEEE)	
Standards Applied:	BS EN 60601-1:2006+A1:2013	Medical Equipment: General requirements	
- -	BS EN 60601-1-2:2015	Medical Equipment: Electromagnetic disturbances	
	BS EN 60601-1-6:2010	Medical Equipment: Usability	
	BS EN 60601-1-11:2015	Medical Equipment: Home healthcare	

We hereby, under our sole responsibility declare, that the product listed above, is in conformity with the Medical Device Regulation (MDR) 2017/745 as a class I medical device accessory based on Annex 8 and the relevant European harmonised standards listed above.

The product concerned has been manufactured under a quality management system (ISO 13485:2016) according to Annex 9 of MDR 2017/745 EU.

This EU declaration of conformity was written accordance to Annex 4 of the MDR, and all supporting documentation is retained at the premises of the manufacturer.

Refer to page 2 for models covered.

Date: 14/12/22

Place: Leicester, UK.

Name: Joe Rider

Position: Director

Version: 1

Q1157 P140 EU Declaration of Conformity

Page 1 of 2



## **EU Declaration of Conformity**

## Models covered by this certificate;

Model:	Description:	Basic UDI-DI (GMN):	UDI-DI (GTIN-14):
P140AA	Movement Sensor: Light Weight, Alert-iT.	506053070P140AR5	50605307027184
P140BA	Movement Sensor: Medium Weight, Alert-iT.	506053070P140BR7	50605307027252
P140CA	Movement Sensor: Heavy Weight, Alert-iT.	506053070P140CR9	50605307027320
P140AB	Movement Sensor: Light Weight, Leisner	506053070P140AR5	50605307027498
P140BB	Movement Sensor: Medium Weight, Leisner.	506053070P140BR7	50605307027566
P140CB	Movement Sensor: Heavy Weight, Leisner.	506053070P140CR9	50605307027634
P140AC	Movement Sensor: Light Weight, Picomed.	506053070P140AR5	50605307027702
P140BC	Movement Sensor: Medium Weight, Picomed.	506053070P140BR7	50605307027870
P140CC	Movement Sensor: Heavy Weight, Picomed.	506053070P140CR9	50605307027948