

Any Nurse or Warden Call system can now become the heart of an advanced Care Monitoring System by the addition of any of the many Alert-iT monitors


Bed or Chair leaving



Incontinence & Vomit



Touch Sensitive Call Buttons



Noise sensitive Call Buttons




General mobility



The radio based care monitoring systems from Alert-iT are recognised as the most advanced and reliable available. They are to be found in care homes and schools as well as many domestic residences. They are specified by local authorities, PCT's, Housing Associations and Assistive Technology providers. The same powerful monitors are all now available for direct connection to established Nurse Call systems, turning the old patient call button facility into a modern, effective care system that improves vigilance for vulnerable people and releases staff for care duties rather than routine checking.

## Radio Gateway



Where direct wired connection is not possible then the new P155C Send-iT Gateway will allow radio based monitors to trigger the Nurse Call system. This will cover our range of fall alarms and call pendants, door sensors, smoke or gas alarms. It can also provide a wirefree and discrete room installation for those residents who may be disturbed by the sight of technology.



Epilepsy & Seizures

